



**MEISINGER**  
IMPLANTS

# SAFESCRAPER TWIST

## Advantages of the disposable semicircular blade

excellent control during the  
harvesting procedure

chamber capacity

straight 3,0 cm<sup>3</sup>  
curved 2,5 cm<sup>3</sup>

maximum cutting efficiency  
with 160° cutting area range

minimally invasive, easy to use  
and well tolerated by the patient

effective on any bone  
surface (plane, concave, convex)

button to prevent accidental  
opening during use or at device disposal



### SAFESCRAPER® TWIST



1 x packaging unit 3



3598  
Safescraper® TWIST  
Straight Version



3987  
Safescraper® TWIST  
Curve Version

