

BONE LEVEL

SURGERY	NC				RC							
	Ø 3.3		ROX-CERA Ø 3.3		Ø 3.75		Ø 4.1		Ø 4.8			
Implant	22009 L 8.0 22010 L 10.0 22011 L 12.0 22012 L 14.0	22200 L 8.0 22201 L 10.0 22202 L 12.0	22050 L 8.0 22051 L 10.0 22052 L 12.0	22201 L 8.0 22202 L 10.0 22203 L 12.0 22204 L 14.0	22060 L 8.0 22061 L 10.0 22062 L 12.0	22005 L 8.0 22006 L 10.0 22007 L 12.0 22008 L 14.0						
Round bur Pilote drill Twist drills												
Twist drills												
Parallel pins Depth gauges												
Depth gauges												
Profile drills												
Taps												
Cover screws	GH 0.0 22013	GH 0.0 22014R*	GH 0.5 22014	GH 0.5 22099	GH 0.5 22099R*	GH 0.0 22098						
Healing Abutments	GH 2.0 Ø 3.6 GH 3.5 22018* GH 5.0 22019	GH 2.0 Ø 4.8 GH 3.5 22021* GH 5.0 22022	GH 2.0 Ø 5.5 GH 4.0 22070 GH 6.0 22071 GH 5.0 22072	GH 2.0 Ø 5.0 GH 4.0 22091* GH 6.0 22092	GH 2.0 Ø 6.5 GH 4.0 22093* GH 6.0 22095	GH 2.0 Ø 7.0 GH 4.0 22123						

BONE LEVEL TAPERED

SURGERY	NC				RC							
	Ø 3.3		Ø 3.75		Ø 4.1		Ø 4.8					
Implants	23001 L 8.0 23002 L 10.0 23003 L 12.0 23004 L 14.0	23005 L 8.0 23006 L 10.0 23007 L 12.0	24001 L 8.0 24002 L 10.0 24003 L 12.0 24004 L 14.0	24005 L 8.0 24006 L 10.0 24007 L 12.0								
Round bur Pilote drill Twist drills												
Twist drills												
Parallel pins Depth gauges												
Depth gauges												
Profile drills												
Taps												
Cover screws	GH 0.0 22013	GH 0.5 22014	GH 0.5 22099	GH 0.0 22098								
Healing Abutments	GH 2.0 Ø 3.6 GH 3.5 22018 GH 5.0 22019	GH 2.0 Ø 4.8 GH 3.5 22020 GH 5.0 22021	GH 2.0 Ø 5.5 GH 4.0 22070 GH 6.0 22071 GH 5.0 22072	GH 2.0 Ø 5.0 GH 4.0 22090 GH 6.0 22092	GH 2.0 Ø 6.5 GH 4.0 22093 GH 6.0 22095	GH 2.0 Ø 7.0 GH 4.0 22123						



Who's who?

OKTAGON®
BONE LEVEL




BONE LEVEL (BL-BLT)

Transfer Parts




NARROW CONNECT (NC)

REGULAR CONNECT (RC)

Impression Caps
Positioning Cylinders





Analogs




Temporary
Abutments





Abutments



Copings



CAD/CAM



Matrix



Additional Parts
Screws

