


TISSUE LEVEL



Who's who?

OKTAGON®
TISSUE LEVEL



CHIRUGIA	RP							WP	
Impianti	 Ø 4.8 Ø 3.3 72186 L 8.0 72187 L 10.0 72199 L 12.0 72267 L 14.0	 Ø 4.8 Ø 3.3 72200 L 8.0 72201 L 10.0 72202 L 12.0	 Ø 4.8 Ø 3.75 31140 L 8.0 31141 L 10.0 31142 L 12.0	 Ø 4.8 Ø 4.1 31149 L 6.0 31151 L 8.0 31152 L 10.0 31153 L 12.0 82647 L 14.0	 Ø 4.8 Ø 4.1 31169 L 8.0 31170 L 10.0 31171 L 12.0	 Ø 4.8 Ø 4.1 38150 L 10.0 38415 L 12.0	 Ø 4.8 Ø 4.1 31147 L 8.0 31150 L 10.0 31175 L 12.0	 Ø 4.8 Ø 4.8 50124 L 6.0 50125 L 8.0 50347 L 10.0 50348 L 12.0	 Ø 6.5 Ø 4.8 93727 L 8.0 93728 L 10.0 93729 L 12.0
Fresa a rosetta Fresa pilota Fresa a spirale									
Fresa a spirale									
Perni di parallelizzazione Calibri di profondità									
Calibri di profondità									
Fresa svasata									
Maschiatori									
Viti di chiusura									
Abutment di guarigione									

Capette per impronta
Cilindri di posizionamento




Analoghi





Abutment provvisori





Abutment



CAD/CAM



Parti aggiuntive
Viti



WIDE PLATFORM (WP)

